AquaShield® system
Système AquaShield®
AquaShield® System
Sistema AquaShield®
Sistema AquaShield®
Sistema AquaShield®
AquaShield® system
AquaShield® systeem
AquaShield® sistemi
AquaShield® 시스템

Reorder No. 00711540
Reorder No. 00711542
Référence de commande 00711540
Référence de commande 00711542
Nachbestell-Nr. 00711540
Nachbestell-Nr. 00711542
N. di riordino 00711540
N. di riordino 00711542
N° de pedido 00711540
N° de pedido 00711542
N.º de encomenda 00711540
N.º de encomenda 00711542
Genbestillingsnr. 00711540
Genbestillingsnr. 00711542
Nabestelnr. 00711540
Nabestelnr. 00711542
Yeni Sipariş Numarası: 00711540
Yeni Sipariş Numarası: 00711542
주문 번호 00711540
주문 번호 00711542

INSTRUCTIONS FOR USE

Instructions d'utilisation
Gebrauchsanweisung
Istruzioni per l'uso
Instrucciones de uso
Instruções de Utilização
Brugsanvisning
Gebruiksaanwijzing
Kullanım talimatları
사용 설명서

*listening...and delivering solutions**

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639.4494

fax +1 440 / 639.4495

email global@usendoscopy.com

www.usendoscopy.com

INTENDED USE:

The AquaShield® system is intended to be used with an air source from an endoscope with the purpose of supplying sterile water to the endoscope during endoscopic procedures. It is compatible with U.S. commercially available sterile water bottles. This device is intended to be used for not more than 24 hours.

For use with Olympus endoscopes series 140, 160, 180, 190, 240, 260 and 290, specifically those endoscopes that have a distal air/water nozzle.

The AquaShield® system (#00711540) is compatible with commercially available water bottles such as Baxter, Hospira and BBraun.

The AquaShield® system (#00711542) is compatible with commercially available water bottles such as Versylene® Fresenius, Baxter Deutschland GmbH, Aqua B. Braun and Versol.

WARNINGS AND PRECAUTIONS:

- This disposable medical device is intended to be used for no more than 24 hours after package has been opened. This device should be discarded within this time period regardless of number of uses.
- Never attach or reattach the AquaShield® system to an endoscope that is/was inside the patient and has not been reprocessed.
- US Endoscopy did not design this device to be reprocessed or used for more than the 24 hour period after the package has been opened, and therefore cannot verify that reprocessing can clean and/or sterilize or maintain the structural integrity of the device to ensure patient and/or user safety.
- Monitor the water level in the water bottle during the procedure. Water feeding will stop once the water level drops below the water bottle tubing.
- Any institution, practitioner, or third party who reprocesses, refurbishes, remanufactures, reesterilizes, and/or reuses this disposable medical device must bear full responsibility for their safety and effectiveness.
- Endoscopic procedures should only be performed by persons having adequate training and familiarity with endoscopic techniques.
- Consult the medical literature relative to complications, hazards, and techniques prior to the performance of any endoscopic procedure.
- Use proper aseptic techniques to avoid contamination of the AquaShield® system during device setup, replacement of water bottles and between patient uses.

CONTRAINDICATIONS:

Contraindications include those specific to any endoscopic procedure.

INSTRUCTIONS FOR USE:

This product is manufactured not to contain latex.

1. Prior to clinical use familiarize yourself with the device.
 - Read the "Warnings and Precautions" and the "Instructions for Use." Review the diagram.
2. Open a bottle of sterile water (250cc, 500cc, or 1000cc).
3. Open the sterile pouch and remove the AquaShield® system. Ensure that the Pinch Clamp is in the open position. Avoid contact with the Water Bottle Tubing to reduce the risk of contamination.
4. Insert the Water Bottle Tubing into the sterile water bottle and firmly tighten the threaded Lid Cap to ensure a seal.
5. Connect the Endoscope Connector to the air/water port connection on the endoscope.
6. Turn on the light source of the processor.
7. **Prime the air/water channel prior to insertion of the endoscope into a patient. If water pressure is low, check to see that the Lid Cap is tightly closed and the Endoscope Connector is firmly attached to the air/water port connection.**
8. **If water bottle needs to be replaced, use proper aseptic technique.**
9. When procedure is completed, close the Pinch Clamp and disconnect the AquaShield® system from the endoscope. This will prevent water from leaking from the Endoscope Connector.
10. Turn off air processor **after** disconnecting the AquaShield® system from the endoscope.
11. Ensure that the device is discarded within 24 hours after the package has been opened. A label is provided with the device.

Product Disposal:

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Issued Date: June 2014

Warning:


An issued or revision date for these instructions is included for the user's information. In the event that two years have elapsed between this date and product use, the user should contact US Endoscopy to determine if additional information is available.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

US Endoscopy, a wholly owned subsidiary of STERIS Corporation

Olympus® is a registered trademark of Olympus Optical Co., Ltd.

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Listening...and delivering solutions®, le logo US Endoscopy listening...and delivering solutions, ainsi que toutes les marques portant le signe ® ou ™ sont déposées auprès de l'U.S. Patent and Trademark Office ou sont des marques commerciales détenues par US Endoscopy.
Listening...and delivering solutions®, das US Endoscopy „listening...and delivering solutions“-Design und alle mit ® oder ™ bezeichneten Marken sind eingetragene Marken, die beim U.S. Patent and Trademark Office (US-Patent- und Markenamt) eingetragen wurden, oder Marken, die sich im Besitz von US Endoscopy befinden.
Listening...and delivering solutions®, il design dello slogan listening...and delivering solutions di US Endoscopy e tutti i marchi contrassegnati da ® o ™ sono stati registrati presso l'US Patent and Trademark Office degli Stati Uniti o sono marchi di fabbrica di proprietà di US Endoscopy.
Listening...and delivering solutions®, o design US Endoscopy listening...and delivering solutions e todas as marcas assinaladas com ® ou ™ estão registadas no Departamento de Registos de Patentes dos EUA (U.S. Patent and Trademark Office) ou são marcas comerciais detidas pela US Endoscopy.
Listening...and delivery solutions®, designet US Endoscopy listening...and delivery solutions og alle varemærker markeret med ® eller ™ er registreret i USA's patent- og varemærkeregister (U.S. Patent and Trademark Office) eller er varemærker, der tilhører US Endoscopy.
Listening...and delivering solutions®, het ontwerp listening...and delivering solutions van US Endoscopy en alle aanduidingen met ® of ™ zijn gedeponneerd bij het Patent and Trademark Office in de Verenigde Staten of zijn handelsmerken van US Endoscopy.