## Biopsy Forceps

### Biopsy Forceps

Pince à biopsie  
Biopsiezange  
Pinza per biopsia  
Fórceps para biopsia  
Pinça de biopsia  
Biopsitang  
Biopsietang  
Biyopsi Forsepsi

Reorder No. 00711201 (230cm)  
Reorder No. 00711202 (230cm)  
Reorder No. 00711203 (160cm)  
Reorder No. 00711204 (160cm)

### Biopsy Forceps – enteroscope

Pince à biopsie - entéroscope  
Biopsiezange - Enteroskop  
Pinza per biopsia - enteroscopio  
Fórceps para biopsia – enteroscopia  
Pinça de biopsia - enteroscópio  
Biopsitang - enteroskop  
Biopsietang - enteroscoop  
Biyopsi Forsepsi - enteroskop

Reorder No. 00711291 (350cm)

## Biopsy Forceps – sheathed

Pince à biopsie - gainée  
Biopsiezange – ummantelt  
Pinza per biopsia - rivestita  
Fórceps para biopsia - con protección  
Pinça de biopsia - revestida  
Biopsitang – sammenfoldet  
Biopsietang – bekleed  
Biyopsi Forsepsi – kılıflı

Reorder No. 00711205 (230cm)

## Biopsy Forceps – hot

### Biopsy Forceps – hot

Pince à biopsie – chaude  
Biopsiezange - für Hochfrequenzstrom  
Pinza per biopsia - a caldo  
Fórceps para biopsia – caliente  
Pinça de biopsia - quente  
Biopsitang - varm  
Hotbiopsietang  
Biyopsi Forsepsi – hot

Reorder No. 00711211 (230cm)

### Biopsy Forceps – hot – enteroscope

Pince à biopsie – chaude - entéroscope  
Biopsiezange - für Hochfrequenzstrom - Enteroskop  
Pinza per biopsia – a caldo - enteroscopio  
Fórceps para biopsia - caliente – enteroscopia  
Pinça de biopsia - quente - enteroscópio  
Biopsitang - varm - enteroskop  
Hotbiopsietang - enteroscoop  
Biyopsi Forsepsi - hot - enteroskop

Reorder No. 00711295 (350cm)

## Biopsy Forceps – pediatric

Pince à biopsie - pédiatrique  
Biopsiezange - Pädiatrie  
Pinza per biopsia - pediatriche  
Fórceps para biopsia - uso pediátrico  
Pinça de biopsia - pediátrica  
Biopsitang - pædiatrisk  
Biopsietang – kind  
Biyopsi Forsepsi - pediyatrik

Reorder No. 00711261 (160cm)  
Reorder No. 00711262 (160cm)

---

## INSTRUCTIONS FOR USE

Instructions d'utilisation  
Gebrauchsanweisung  
Istruzioni per l'uso  
Instrucciones de uso  
Instruções de Utilização  
Brugsanvisning  
Gebruiksaanwijzing  
Kullanım talimatları

---

*listening...and delivering solutions\**

5976 Heisley Road  
Mentor, OH 44060 USA

phone +1 440 / 639.4494

fax +1 440 / 639.4495

email [global@usendoscopy.com](mailto:global@usendoscopy.com)

[www.usendoscopy.com](http://www.usendoscopy.com)

Rx Only (U.S.A.)

**INTENDED USE:**

Disposable biopsy forceps are used to obtain mucosal tissue samples during endoscopic procedures.

**WARNINGS AND PRECAUTIONS:**

- Endoscopic biopsy should only be performed by persons having adequate training and familiarity with endoscopic biopsy techniques. Consult the medical literature relative to techniques, contraindications, complications and hazards prior to any endoscopic procedure. If reference information for endoscopic technique or contraindications is needed call US Endoscopy.
- Do not use this device if there is any apparent shipping or handling damage.
- Endoscopic biopsy should not be attempted until proficiency has been developed with biopsy forceps.
- Excessive force should not be necessary and should be avoided as not to damage the forceps.
- Biopsy forceps should never be forced through the endoscope's channel. If resistance is met due to excess angulation of the endoscope, it may be necessary to decrease the angulation to allow passage of the forceps.
- Biopsy forceps should not be extracted through the biopsy channel unless the biopsy cups have been completely closed (with or without a specimen.)
- US Endoscopy did not design this device to be reprocessed or reused, and therefore cannot verify that reprocessing can clean and/or sterilize or maintain the structural integrity of the device to ensure patient and/or user safety.
- **These disposable medical devices are not intended for reuse.** Any institution, practitioner or third party who reprocesses, refurbishes, remanufactures, sterilizes and/or reuses these disposable medical devices must bear full responsibility for product safety and effectiveness.
- **For product numbers 00711211, 00711295: Endoscopic hot biopsy** should only be performed by persons having a thorough understanding of the principals of diathermic energy. Consult the electrosurgical generator manufacturer's instruction booklet for proper settings and use of the generator. In order to ensure that the insulating properties of the device are not compromised, do not exceed the maximum rated peak recurring voltage stated for each high frequency mode (see chart below).

HIGH FREQUENCY MODE	MAXIMUM RATED RECURRING PEAK VOLTAGE
Cut	1200 Volts
Coagulation Burst	2100 Volts
Coagulation Spray	5200 Volts

- Care should be exercised when using hot biopsy forceps to avoid causing unwanted burning or perforation of the wall of the GI tract.

**DIRECTIONS FOR USE:**

1. Select the forceps compatible with the biopsy channel of the endoscope.
2. Open the sterile pouch. Remove and inspect forceps for obvious damage. Uncoil the forceps and actuate the handle several times to be sure the biopsy cups open and close properly. Should there be evidence of damage, or if the unit is not functioning properly, **DO NOT USE THIS PRODUCT.**
3. After endoscopic visualization of the area for biopsy, place thumb in loop of forceps handle, rest middle and index finger (placing gentle pressure) on the handle's spool edge.
4. With gentle pressure on the handle to keep the cups closed, insert the forceps through the biopsy channel of the endoscope using short strokes until the cups are visualized exiting the endoscope.
5. Maneuver the forceps towards the targeted site. Open the cups by activating the handle. Advance the opened forceps over the targeted site and close the forceps. Use only enough pressure to bite the tissue. Over-exertion could cause forceps to become misaligned or to fail.
6. Continue to apply gentle pressure on the handle as the forceps are withdrawn from the channel.
7. Remove and prepare the tissue specimen according to standard technique for histologic evaluation.

**DIRECTIONS FOR USE: (FOR PRODUCT NUMBERS 00711211, 00711295 Only)**

1. Follow steps 1 – 4 above.
2. Attach the correct active cord to the handle and plug into generator. Set the generator to the pre-determined settings.
3. Maneuver the forceps towards the targeted site. Open the cups by activating the handle. Advance the opened forceps over the targeted site and close the forceps. Use only enough pressure to bite the tissue. Over-exertion could cause forceps to become misaligned or to fail.
4. The time and power settings should be determined clinically to control hemostasis.
5. Once the biopsy has been completed and hemostasis has been attained, the procedure can continue to the next step.
6. Continue to apply gentle pressure on the handle as the forceps are withdrawn from the channel.
7. Remove and prepare the tissue specimen according to standard technique for histologic evaluation.

**PRODUCT DISPOSAL:**

After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

**ISSUED DATE:** April 2010

**Warning:**

An issued or revision date for these instructions is included for the user's information. In the event two years has elapsed between this date and product use, the user should contact US Endoscopy to determine if additional information is available.

Listening...and delivering solutions<sup>®</sup>, the US Endoscopy listening...and delivering solutions design, and all marks denoted with <sup>®</sup> or <sup>™</sup> are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy

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**Listening...and delivering solutions®**, le logo US Endoscopy listening...and delivering solutions, ainsi que toutes les marques portant le signe ® ou ™ sont déposées auprès de l'U.S. Patent and Trademark Office ou sont des marques de commerce détenues par US Endoscopy.
**Listening...and delivering solutions®**, das US Endoscopy „listening...and delivering solutions“-Design und alle mit ® oder ™ bezeichneten Marken sind eingetragene Marken, die beim U.S. Patent and Trademark Office (US-Patent- und Markenamt) eingetragen wurden, oder Marken, die sich im Besitz von US Endoscopy befinden.
**Listening...and delivering solutions®**, o design US Endoscopy listening...and delivering solutions e todas as marcas assinaladas com ® ou ™ estão registadas no Departamento de Registos de Patentes dos EUA (U.S. Patent and Trademark Office) ou são marcas comerciais detidas pela US Endoscopy
**Listening...and delivering solutions®**, US Endoscopy listening...and delivering solutions design og alle varemærker markeret med ® eller ™ er registreret i USA's patentregister (U.S. Patent and Trademark Office) eller er varemærker, der tilhører US Endoscopy.
**Listening...and delivering solutions®**, het ontwerp listening...and delivering solutions van US Endoscopy en alle aanduidingen met ® of ™ zijn geregistreerd bij het Patent and Trademark Office in de Verenigde Staten of zijn handelsmerken van US Endoscopy.
**Listening...and delivering solutions®**, US Endoscopy listening...and delivering solutions의 디자인과 • 또는 ™ 이 표시된 모든 마크는 미국 특허청(the U.S. Patent and Trademark Office)에 등록되었거나 US Endoscopy 소유의 등록 상표입니다.