

# Advanced Tissue Resection

iSnare™ Case Report Series • Report 3



**“The iSnare is an excellent device. It saved the patient from a large surgery. I liked the stiffness of the snare and time savings of 5 to 10 minutes in the procedure.”**

*Jeffrey L. Ponsky, M.D.  
Oliver H. Payne Professor  
and Chairman,  
Department of Surgery,  
CWRU and University  
Hospitals of Cleveland*

## Procedure

- Colonoscopy with resection of large sessile cecal polyp

## Indications

- Colonic polyps

## Findings

The colon was normal except for a previously noted 4cm sessile polyp just distal to the ileo-cecal valve on the lateral wall.

## Endoscopic Mucosal Resection


The previously noted 4cm sessile polyp was identified. Using the iSnare™, a multi-functional injection and snaring device, the lesion was injected submucosally with saline a total of 5 times. The device's snare was used to resect the lifted polyp in a piecemeal fashion. There was no bleb loss due to absorption because instrument exchange was not necessary. The entire lesion was resected. Using the Roth Net® retrieval device, most of the specimen fragments were retrieved for pathological examination.



1. Submucosal bleb formation using normal saline



2. Snare deployment prior to transection



3. Polyp fragment retrieval post resection



*listening...and delivering solutions*

5976 Heisley Road  
Mentor OH 44060

phone 440 / 639,4494

fax 440 / 639,4495

customer service 800 / 769,8226

[www.usendoscopy.com](http://www.usendoscopy.com)

This case report is not intended to provide device or clinical instruction or intended to be a substitute for the Instructions For Use (IFU) which are packaged with the device. Please consult the IFU prior to using the iSnare™ device.

iSnare™ is a trademark of US Endoscopy, patents pending.

Roth Net® is a registered U.S. trademark of US Endoscopy, U.S. patent numbers 6,814,739 and 5,643,283, other patents pending.

760506 Rev. A