eTrap[®] – polyp trap

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eTrap[®] – Polypenfalle

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Rx Only (U.S.A.)

Intended Use:

The single use eTrap[®] – polyp trap is used for suction retrieval of endoscopically removed polyps.

Warnings and Precautions:

- US endoscopy did not design this device to be reprocessed or reused, and therefore cannot verify that reprocessing can clean and/or sterilize or maintain the structural integrity of the device to ensure patient and/or user safety.
- This disposable medical device is not intended for reuse. Any institution, practitioner, or third party who reprocesses, refurbishes, remanufactures, reesterilizes, and or reuses this disposable medical device must bear full responsibility for its safety and effectiveness.
- Two specimen retrieval strainers accompany the eTrap[®]. Always insert the extra strainer into the eTrap[®] upon removal of the first strainer; this practice will guard against inadvertent suctioning of a polyp through the eTrap[®] and into the main suction system.

Directions for Use:

1. Remove both specimen retrieval strainers from the packaging and insert one of the strainers into the rectangular side port of the eTrap[®]. Keep the second specimen retrieval strainer available for use as described in steps 6 and 7.
2. Connect the tube from the top of the eTrap[®] to the endoscope's suction barb.
3. Connect the main suction canister's tubing to the eTrap[®] device's outgoing suction barb.
4. Suction through the endoscope using standard endoscopic technique until the excised polyp is entrapped in the specimen retrieval strainer. The magnifying window on top of the eTrap[®] aids in visualization.
5. Remove the specimen retrieval strainer from the eTrap[®].
6. Insert the extra specimen retrieval strainer into the eTrap[®] for additional usage.
7. The process of removing and replacing the specimen retrieval strainers in tandem (as needed), prevents loss of suction and allows the availability to retrieve additional polyps while the previously retrieved polyp is prepared by clinical staff. Preparation of the retrieved polyp is achieved by removing it from the specimen retrieval strainer and placing it in the appropriate fixative for histopathology in accordance with institutional guidelines.
8. When the eTrap[®] is no longer needed, disconnect the main suction canister's tubing from the eTrap[®], disconnect the eTrap[®] device's tubing from the endoscope and reconnect the main suction canister's tubing to the endoscope.

Product Disposal:

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Packaging is constructed of polyethylene terephthalate. Please recycle.

Issued Date: April 2010

Warning:

An issued or revision date for these instructions is included for the user's information. In the event that two years have elapsed between this date and product use, the user should contact US Endoscopy to determine if additional information is available.

Listening...and delivering solutions[®], the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

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Listening...and delivering solutions®, le logo US Endoscopy listening...and delivering solutions, ainsi que toutes les marques portant le signe ® ou ™ sont déposées auprès de l'U.S. Patent and Trademark Office ou sont des marques commerciales détenues par US Endoscopy.
Listening...and delivering solutions®, das US Endoscopy „listening...and delivering solutions“-Design und alle mit ® oder ™ bezeichneten Marken sind eingetragene Marken, die beim U.S. Patent and Trademark Office (US-Patent- und Markenamt) registriert wurden, oder Marken, die sich im Besitz von US Endoscopy befinden.
Listening...and delivering solutions®, o design US Endoscopy listening...and delivering solutions e todas as marcas assinaladas com ® ou ™ estão registadas no Departamento de Registos de Patentes dos EUA (U.S. Patent and Trademark Office) ou são marcas comerciais detidas pela US Endoscopy
Listening...and delivering solutions®, US Endoscopy listening...and delivering solutions design og alle varemærker markeret med ® eller ™ er registreret i U.S. Patent and Trademark Office eller er varemærker, der tilhører US Endoscopy
Listening...and delivering solutions[®], het ontwerp listening...and delivering solutions van US Endoscopy en alle aanduidingen met [®] of [™] zijn geregistreerd bij het Patent and Trademark Office in de Verenigde Staten of zijn handelsmerken van US Endoscopy.
Listening...and delivering solutions®, US Endoscopy listening...and delivering solutions 디자인, 그리고 ® 또는 ™으로 표시된 모든 상표는 미국 특허청(U.S. Patent and Trademark Office)에 등록되어 있거나 US Endoscopy가 소유하고 있는 상표입니다.