Guardus® overtube – esophageal – disposable

Overtube Guardus® - oesophagien - usage unique
Guardus® Overtube – Ösophagus – Einweg
Guardus® - overtube - esofageo – monouso
Guía esofagal Guardus® – desechable
Overtube Guardus® – esofágico – descartável
Guardus® overtube – øsofagus – engangsartikel
Guardus® overtube – oesofageaal – wegwerpbaar
Guardus® overtube – özofajeal- tek kullanımlık
Guardus® 오버튜브 - 식도 - 일회용

Reorder No: 00711146 / SI711146 / 00711147 / SI711147

Référence de commande 00711146 / SI0071146 / 00711147 / SI711147
Nachbestell-Nr. 00711146 / SI0071146 / 00711147 / SI711147
N. di riordino 00711146 / SI0071146 / 00711147 / SI711147
Nº de pedido 00711146 / SI0071146 / 00711147 / SI711147
N.º de encomenda 00711146 / SI0071146 / 00711147 / SI711147
Genbestillingsnr. 00711146 / SI0071146 / 00711147 / SI711147
Nabestelnr.00711146 / SI0071146 / 00711147 / SI711147
Yeni Sipariş Numarası: 00711146 / SI0071146 / 00711147 / SI711147
주문 번호 00711146 / SI0071146 / 00711147 / SI711147

Guardus® overtube – gastric – disposable

Overtube Guardus® - gastrique - usage unique
Guardus® Overtube – Magen – Einweg
Guardus® - overtube - gastrico – monouso
Guía gástrica Guardus® – desechable
Overtube Guardus® – gástrico – descartável
Guardus® overtube – gastrisk – engangsartikel
Guardus® overtube – voor maag – wegwerpbaar
Guardus® overtube – gastrik- tek kullanımlık
Guardus® 오버튜브 - 위장 - 일회용

Reorder No. 00711148 / SI711148 / 00711149 / SI711149

Référence de commande 00711148 / SI711148 / 00711149 / SI711149
Nachbestell-Nr. 00711148 / SI711148 / 00711149 / SI711149
N. di riordino 00711148 / SI711148 / 00711149 / SI711149
Nº de pedido 00711148 / SI711148 / 00711149 / SI711149
N.º de encomenda 00711148 / SI711148 / 00711149 / SI711149
Genbestillingsnr. 00711148 / SI711148 / 00711149 / SI711149
Nabestelnr. 00711148 / SI711148 / 00711149 / SI711149
Yeni Sipariş Numarası: 00711148 / SI711148 / 00711149 / SI711149
주문 번호 00711148 / SI711148 / 00711149 / SI711149

INSTRUCTIONS FOR USE

Instructions d'utilisation
Gebrauchsanweisung
Istruzioni per l'uso
Instrucciones de uso
Instruções de Utilização
Brugsanvisning
Gebruiksaanwijzing
Kullanım talimatları
사용 설명서

listening... and delivering solutions®

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639.4494

fax +1 440 / 639.4495

email global@usendoscopy.com

www.usendoscopy.com

NON-STERILE

Rx Only (U.S.A.)

Reorder No. 00711146 / SI711146

(25cm length, 16.7mm ID, 19.5mm OD, endoscope OD range 8.6mm-10.0mm)

Reorder No. 00711147 / SI711147

(25cm length, 16.7mm ID, 19.5mm OD, endoscope OD range 10.0mm-11.7mm)

Reorder No. 00711148 / SI711148

(50cm length, 16.7mm ID, 19.5mm OD, endoscope OD range 8.6mm-10.0mm)

Reorder No. 00711149 / SI711149

(50cm length, 16.7mm ID, 19.5mm OD, endoscope OD range 10.0mm-11.7mm)

Overtube assembly contains (1) inner tube with tapered tip and (1) overtube with tethered insufflation cap

Intended Use:

- The disposable overtube is a device used in conjunction with a flexible endoscope for foreign body or tissue retrieval, and/or for endoscopic procedures requiring multiple endoscope intubations.

Warnings and Precautions:

- Endoscopic procedures should only be performed by persons having adequate training and familiarity with endoscopic techniques.
- Consult the medical literature relative to complications, hazards and techniques prior to the performance of any endoscopic procedure.
- Use of an overtube, as well as foreign body retrieval, should not be attempted unless proficiency in technique has been developed by the clinician.
- Baseline EGD is recommended prior to overtube usage to minimize potential patient complications.
- Use of overtubes has been associated with esophageal or pharyngeal laceration or perforation. Such complications should be monitored for during and after the procedure and treated appropriately.
- Ensure that the OD of the endoscope is compatible with the ID of the inner tube with tapered tip (see product number size specifications).
- Performance may be compromised if proper lubrication is not applied to both the inner and outer surfaces of the overtube components and endoscope.
- Once the overtube is in place, it should not be advanced beyond the endoscope. Any migration or advancement could result in tissue entrapment and mucosal injury. If repositioning is required, please follow directions as listed in steps 1-4 in the "Repositioning the Overtube" section.
- US Endoscopy did not design this device to be reprocessed or reused, and therefore cannot verify that reprocessing can clean and/or sterilize or maintain the structural integrity of the device to ensure patient and/or user safety.
- This disposable medical device is not intended for reuse. Any institution, practitioner, or third party who reprocesses, refurbishes, remanufactures, reesterilizes, and/or reuses this disposable medical device must bear full responsibility for its safety and effectiveness.

Band Ligation Devices:

- The overtube's insufflation cap **must be removed** to facilitate passage of band ligation devices. **The clinician must ensure size compatibility between specific manufacturers' band ligation device diameters and the ID (inner diameter) of the overtube prior to concomitant use.** Be aware that the overtube must be successfully placed according to these Instructions for Use prior to the use of a band ligation device.

Contraindications:

- Esophageal bleeding and/or laceration; laryngeal perforation; trauma to teeth, gums and/or pharynx; aspiration pneumonia.
- Contraindications also include those specific to any endoscopic procedure, use of an overtube and/or foreign body retrieval.

Pre-Use Instructions:

Prior to clinical use, familiarize yourself with the device.

- Read the "Instructions for Use" and review the photographs.
- Inspect the package for shipping or handling damage. If damage is evident, do not use this device.
- For user convenience, the inner tube handles have been color-coded. The yellow inner tube fits endoscope OD (outer diameter) range 8.6mm-10.0mm and the orange inner tube fits endoscope OD (outer diameter) range 10.0-11.7mm. Ensure that the OD (outer diameter) of the endoscope is compatible with the ID (inner diameter) of the inner tube (see product number color and size specifications) prior to opening package.
- This device is packaged in two pieces – the *inner tube with tapered tip*, and the *overtube with tethered insufflation cap*. Prior to use, remove lubrication tag from inner tube and overtube. **Generously lubricate both components with a water soluble lubricant** as described below in the Directions for Use. Do not assemble without generous lubrication.

Directions for Use:

- Prepare patient for endoscopy. It is recommended that a 60FR bite block be used to prevent inadvertent collapse of the overtube by patient gums or dentition.
- After removing the overtube assembly components from the packaging, **generously lubricate the inner and outer surfaces of both tubes with a water soluble lubricant.**
- Load the **well lubricated** inner tube *into* the **well lubricated** overtube. (See Fig. 1)
- Lock the inner tube into the overtube by engaging the cap threads into the handle and *twisting the cap clockwise*. (See Fig. 2)
- Once locked, visually inspect the distal insertion end of the engaged dual device to ensure that the tapered tip forms a *gap-free interface* with the ID of the overtube.
- Generously lubricate** the lower one third of the distal (insertion end) of the endoscope, and *backload the well lubricated overtube assembly onto the well lubricated endoscope*, sliding it to the proximal end of the endoscope's insertion tube. (See Fig. 3)
- Ensure that the tapered tip forms a *gap-free interface* with the OD of the endoscope, and that it slides freely over the endoscope. (See Fig. 4)
- Perform baseline EGD.
- If the clinician determines that passage of the overtube assembly is achievable/desired, the overtube assembly should be *gently passed forward* over the **well lubricated** endoscope insertion tube until it reaches the desired anatomical location. Visually observe the overtube assembly depth of insertion by noting the centimeter markings on the overtube.

10. Once the overtube assembly is in the desired location, turn the proximal lock on the inner tube cap *counter-clockwise* to unlock it, **raise the patient's chin**, and *withdraw both the inner tube and endoscope simultaneously* from the patient. You must maintain position of the overtube as you withdraw the inner tube and endoscope. (See Fig. 5)
11. The *tethered insufflation cap* should then be locked onto the overtube by twisting it *clockwise*. (See Fig. 6)
12. Once the endoscope is removed from the inner tube, the endoscopic procedure may be performed through the overtube. (See Fig. 7)

Repositioning the Overtube:

1. If repositioning of the overtube deeper into the gastrointestinal lumen is required, extubate the endoscope from the overtube, remove the insufflation cap, **generously lubricate the inner and outer surfaces of the inner tube**, and *backload* the inner tube with tapered tip onto the **well lubricated** endoscope.
2. Extend the endoscope tip approximately 15 cm (6 inches) past the distal end of the inner tube during repeat intubation to allow for direct endoscopic visualization of the mucosal lining at the distal end of the overtube. This will allow the clinician to endoscopically clear the area (through insufflation) of any enveloped mucosa prior to advancing the inner tube.
3. Once the clinician is sure that the ID of the overtube is free of mucosa, the inner tube may be fully advanced, engaging with the overtube to create a gap-free interface. The clinician can then reposition the overtube assembly to the desired position.
4. Once the desired overtube assembly position is achieved, the clinician should repeat steps 10, 11 and 12.

Product Disposal:

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Issued Date: December 2011

Warning:

An issued or revision date for these instructions is included for the user's information. In the event that two years have elapsed between this date and product use, the user should contact US Endoscopy to determine if additional information is available.

Packaging is constructed of polyethylene terephthalate. Please recycle.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

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Listening...and delivering solutions®, le logo US Endoscopy listening...and delivering solutions, ainsi que toutes les marques portant le signe ® ou ™ sont déposées auprès de l'U.S. Patent and Trademark Office ou sont des marques commerciales détenues par US Endoscopy.
Listening...and delivering solutions®, das US Endoscopy „listening...and delivering solutions“-Design und alle mit ® oder ™ bezeichneten Marken sind eingetragene Marken, die beim U.S. Patent and Trademark Office (US-Patent- und Markenamt) registriert wurden, oder Marken, die sich im Besitz von US Endoscopy befinden.
Listening...and delivering solutions®, o design US Endoscopy listening...and delivering solutions e todas as marcas assinaladas com ® ou ™ estão registadas no Departamento de Registos de Patentes dos EUA (U.S. Patent and Trademark Office) ou são marcas comerciais detidas pela US Endoscopy
Listening...and delivering solutions[®], US Endoscopy listening...and delivering solutions design og alle varemærker markeret med © eller [™] er registreret i U.S. Patent and Trademark Office eller er varemærker, der tilhører US Endoscopy.
Listening...and delivering solutions®, het ontwerp listening...and delivering solutions van US Endoscopy en alle aanduidingen met © of ™ zijn geregistreerd bij het Patent and Trademark Office in de Verenigde Staten of zijn handelsmerken van US Endoscopy.
Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions 디자인, 그리고 ® 또는 ™로 표시된 모든 상표는 미국 특허청 (U.S. Patent and Trademark Office) 에 등록되어 있거나 US Endoscopy가 소유하고 있는 상표입니다.