## Roth Net® – foreign body – enteroscope

Roth Net® – corps étranger - entéroscope  
Roth Net® – Fremdkörper – Enteroskop  
Roth Net® – retino corpi estranei - enteroscopia  
Roth Net® – extracción de cuerpos extraños – enteroscopia  
Roth Net® – corpo estranho – enteroscópio  
Roth Net® – fremmedlegeme – enteroskop  
Roth Net® – vreemd object – enteroscoop  
Roth Net® – yabancı cisim – enteroskop  
Roth Net® – 이물 – 장경

## Roth Net® – foreign body – maxi

Roth Net® – corps étranger - maxi  
Roth Net® – Fremdkörper – Maxi  
Roth Net® – retino corpi estranei - maxi  
Roth Net® – extracción de cuerpos extraños – maxi  
Roth Net® – corpo estranho – maxi  
Roth Net® – fremmedlegeme – maxi  
Roth Net® – vreemd object – maxi  
Roth Net® – yabancı cisim – maksı  
Roth Net® – 이물 - 맥시

## Roth Net® – foreign body – standard

Roth Net® – corps étranger – standard  
Roth Net® – Fremdkörper – Standard  
Roth Net® – retino corpi estranei - standard  
Roth Net® – extracción de cuerpos extraños – estándar  
Roth Net® – corpo estranho – standard  
Roth Net® – fremmedlegeme – standard  
Roth Net® – vreemd object – standaard  
Roth Net® – yabancı cisim – standart  
Roth Net® – 이물 - 표준

## Roth Net® – polyp retrieval

Roth Net® – extraction de polype  
Roth Net® – Polypenentfernung  
Roth Net® – retino per polipi  
Roth Net® – extirpación de pólipos  
Roth Net® – recolha de pólipos  
Roth Net® – polypslynge  
Roth Net® – ophalen poliep  
Roth Net® – polip çıkarma  
Roth Net® - 용종 절제

## Reorder No. 00711152

Référence de commande 00711152  
Nachbestell-Nr. 00711152  
N. di riordino 00711152  
N° de pedido 00711152  
N.º de encomenda 00711152  
Genbestillingsnr. 00711152  
Nabestelnr. 00711152  
Yeni Sipariş Numarası: 00711152  
주문 번호00711152

## Reorder No. 00711059

Référence de commande 00711059  
Nachbestell-Nr. 00711059  
N. di riordino 00711059  
N° de pedido 00711059  
N.º de encomenda 00711059  
Genbestillingsnr. 00711059  
Nabestelnr. 00711059  
Yeni Sipariş Numarası: 00711059  
주문 번호00711059

## Reorder No. 00711050 / 00711052 / 00711053

Référence de commande 00711050 / 00711052 / 00711053  
Nachbestell-Nr. 00711050 / 00711052 / 00711053  
N. di riordino 00711050 / 00711052 / 00711053  
N° de pedido 00711050 / 00711052 / 00711053  
N.º de encomenda 00711050 / 00711052 / 00711053  
Genbestillingsnr. 00711050 / 00711052 / 00711053  
Nabestelnr. 00711050 / 00711052 / 00711053  
Yeni Sipariş Numarası: 00711050 / 00711052 / 00711053  
주문 번호00711050 / 00711052 / 00711053

## Reorder No. 00711150 / 00711151

Référence de commande 00711150 / 00711151  
Nachbestell-Nr. 00711150 / 00711151  
N. di riordino 00711150 / 00711151  
N° de pedido 00711150 / 00711151  
N.º de encomenda 00711150 / 00711151  
Genbestillingsnr. 00711150 / 00711151  
Nabestelnr. 00711150 / 00711151  
Yeni Sipariş Numarası: 00711150 / 00711151  
주문 번호00711150 / 00711151

## INSTRUCTIONS FOR USE

Instructions d'utilisation  
Gebrauchsanweisung  
Istruzioni per l'uso  
Instrucciones de uso  
Instruções de Utilização  
Brugsanvisning  
Gebruiksaanwijzing  
Kullanım talimatları  
사용 설명서

**NON-STERILE**

Rx Only (U.S.A.)

#### Intended Use:

The disposable Roth Net® retrieval device product line is used in the endoscopic retrieval of foreign body, food bolus and excised tissue such as polyps.

<u>Description</u>	<u>Product Number</u>	<u>Sheath Diameter</u>	<u>Sheath Length</u>
Roth Net® – foreign body – standard	00711053	2.5mm	160cm
Roth Net® – foreign body – standard	00711050 / 00711052	2.5mm	230cm
Roth Net® – foreign body – maxi	00711059	3.0mm	160cm
Roth Net® – polyp retrieval	00711150	2.5mm	230cm
Roth Net® – polyp retrieval	00711151	2.5mm	230cm
Roth Net® – foreign body – enteroscope	00711152	2.5mm	350cm

#### Warnings and Precautions:

- Do not use this device if there is any apparent shipping or handling damage.
- The endoscopic retrieval of polyps or foreign bodies should only be performed by persons having adequate training and familiarity with endoscopic techniques. Consult the medical literature relative to techniques, complications and hazards prior to the performance of any endoscopic procedure. If reference information for endoscopic techniques utilizing the Roth Net® retrieval device product line is desired/needed, contact US Endoscopy; case reports describing techniques are available.
- Avoid blindly passing the Roth Net® device past any foreign object, particularly if the entire lumen is blocked.
- Care should be exercised when grasping the targeted polyp/object to be retrieved to avoid inadvertently grasping tissue or organs not intended for retrieval.
- Keep gentle traction on the device during retrieval so that the retrieved object does not become loosened, dislodged or aspirated into the trachea.
- These devices are not recommended for the retrieval of sharp foreign objects.
- The following conditions may not allow the Roth Net® device to function properly: (1) advancing the handle to the open position with too much speed or force, (2) attempting to pass or open the device in an extremely articulated endoscope, (3) attempting to actuate the device in an extremely coiled position and/or (4) actuating the device when the handle is at an acute angle in relation to the sheath.
- Short strokes, 1"-1.5" (2.5cm – 3.8cm) in length, are recommended throughout device passage to avoid sheath kinking.
- Insertion forces may be greater for product #00711152 during enteroscopic use due to scope configuration and scope/device length. The user may experience moderate resistance in device passage as the device exits the distal end of the enteroscope.
- US Endoscopy did not design this device to be reprocessed or reused, and therefore cannot verify that reprocessing can clean and/or sterilize or maintain the structural integrity of the device to ensure patient and/or user safety.
- The Roth Net® retrieval device product line is intended for single patient use only. Any institution, practitioner, or third party who reprocesses, refurbishes, remanufactures, reesterilizes, and/or reuses these disposable devices must bear full responsibility for their safety and effectiveness.

#### Directions for Use:

1. Prior to clinical use you should familiarize yourself with the device.
  - Read the "Warnings and Precautions" and the "Directions for Use." Review the diagrams.
2. **To repack the Roth Net® device into the sheath for introduction into the endoscope:**
  - Uncoil the entire device and drape in a "U" shaped configuration, holding the proximal device end in one hand, and the distal sheath end in the opposite hand.
  - With the device in this position, carefully retract the handle until the Roth Net® device is completely withdrawn into its sheath (see Fig. 3).
3. Check that the accessory channel diameter of the endoscope is compatible with the diameter of the Roth Net® sheath.
4. Once the polyp/object has been endoscopically identified, advance the retracted device into the accessory channel of the endoscope using short strokes until the distal end of the sheath is endoscopically visualized.
5. Advance the distal end of the sheath slightly past the polyp/object to be retrieved, however, avoid blindly passing the device past any polyp/foreign object if the entire lumen is blocked. Open the Roth Net® device by advancing the handle forward until it stops (see Fig. 1, 4). Confirm the net is fully open via endoscopic observation. **The following conditions may cause the device to function improperly: (1) advancing the handle to the open position with too much speed or force, (2) attempting to pass or open the device in an extremely articulated endoscope, (3) attempting to actuate the device in an extremely coiled position and/or (4) actuating the device when the handle is at an acute angle in relation to the sheath.**
6. Endoscopically manipulate the Roth Net® device over the polyp/object. **For best results, the object should be centered and proximal in the net prior to closure.** Once entrapment is accomplished, retract the handle until it stops (see Fig. 2). Continuous gentle traction should be applied on the handle to keep the Roth Net® device closed. Excess pressure may cause the net to rupture or tear.
7. Retract the sheath, retrograde, toward the endoscope so the Roth Net® device is endoscopically observed yet not obscuring the endoscopic view. Endoscopic observation is necessary during extubation as not to lose sight of the polyp/object during removal.
8. Once the endoscope and the Roth Net® have been extubated, the retrieved polyp/object can be removed from the net by advancing the handle forward to open the Roth Net® device (see Fig. 1, 4). Specimens should be prepared according to standard technique for histologic evaluation.
9. **To repack the Roth Net® into the sheath for introduction to, or removal from the endoscope, or for additional retrievals:**
  - Retract the handle until the net of the Roth Net® device is completely retracted into the sheath (see Fig. 3).
10. Once the endoscopic retrieval has been completed, remove the device from the endoscope.

#### Product Disposal:

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

**Issued Date: April 2010**

#### Warning:

An issued or revision date for these instructions is included for the user's information. In the event that two years have elapsed between this date and product use, the user should contact US Endoscopy to determine if additional information is available.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, are trademarks owned by US Endoscopy. U.S. Patent numbers 6,814,739; 5,643,283.

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**Listening...and delivering solutions®**, le logo US Endoscopy listening...and delivering solutions, ainsi que toutes les marques portant le signe ® ou ™ sont déposées auprès de l'U.S. Patent and Trademark Office ou sont des marques commerciales détenues par US Endoscopy. Brevets 6 814 739, 5 643 283, et brevets en attente d'homologation.