Bite Block

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Bite Block – pediatric – latex free strap

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Bite Block – intermediate – latex free strap **Reorder No. 00712809**

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Rx Only (U.S.A.)

Intended Use:

The bite block is used to protect the endoscope insertion tube from being bitten by the patient. Bite blocks are available with and without retaining straps.

Warnings and Precautions:

- When placing the bite block in the patient's mouth, care should be exercised to avoid placing fingers between the bite block and the patient's teeth.
- Flexible endoscopy procedures should only be performed by persons having adequate training and familiarity with endoscopic techniques. Consult the medical literature relative to techniques, complications and hazards prior to the performance of an endoscopic procedure.
- US endoscopy did not design this device to be reprocessed or reused, and therefore cannot verify that reprocessing can clean and/or sterilize or maintain the structural integrity of the device to ensure patient and/or user safety.
- The bite block is intended for single patient use only. Any institution, practitioner, or third party who reprocesses, refurbishes, remanufactures, reesterilizes, and/or reuses this disposable device must bear full responsibility for its safety and effectiveness.

Directions for use:

- The bite block is provided non-sterile; should sterilization be required, cold soak or ethylene oxide (ETO) methods are compatible with polymer. **Bite blocks cannot be flash or steam autoclaved.**
- Place the bite block in the patient's mouth. If the retaining strap is used, secure the strap on the bite block.
- To secure retaining strap (if applicable): Insert the retaining strap through posterior aspect of side port and thread over cleat. If strap is Velcro, press strips together to fasten and pull apart to unfasten.

Product Disposal:

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Issued Date: April 2010

Warning:

An issued or revision date for these instructions is included for the user's information. In the event two years has elapsed between this date and product use, the user should contact US Endoscopy to determine if additional information is available.

Listening...and delivering solutions[®], the US Endoscopy listening...and delivering solutions design, and all marks denoted with [®] or [™] are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy

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Listening...and delivering solutions[®], le logo US Endoscopy listening...and delivering solutions, ainsi que toutes les marques portant le signe[®] ou[™] sont déposées auprès de l'U.S. Patent and Trademark Office ou sont des marques commerciales détenues par US Endoscopy.
Listening...and delivering solutions[®], das US Endoscopy „listening...and delivering solutions“-Design und alle mit [®] oder [™] bezeichneten Marken sind eingetragene Marken, die beim U.S. Patent und Trademark Office (US-Patent- und Markenamt) registriert wurden, oder Marken, die sich im Besitz von US Endoscopy befinden.
Listening...and delivering solutions[®], o design US Endoscopy listening...and delivering solutions e todas as marcas assinaladas com [®] ou [™] estão registadas no Departamento de Registos de Patentes dos EUA (U.S. Patent and Trademark Office) ou são marcas comerciais detidas pela US Endoscopy
Listening...and delivering solutions[®], US Endoscopy listening...and delivering solutions design og alle varemærker markeret med [®] eller [™] er registreret i U.S. Patent and Trademark Office eller er varemærker, der tilhører US Endoscopy
Listening...and delivering solutions[®], het ontwerp listening...and delivering solutions van US Endoscopy en alle aanduidingen met [®] of [™] zijn geregistreerd bij het Patent and Trademark Office in de Verenigde Staten of zijn handelsmerken van US Endoscopy
Listening...and delivering solutions®. US Endoscopy listening...and delivering solutions 디자인, 그리고 ® 또는 ™으로 표시된 모든 상표는 미국 특허청 (U.S. Patent and Trademark Office) 에 등록되어 있거나 US Endoscopy 가 소유하고 있는 상표입니다.