



















US Endoscopy | Roth Net® Retriever | Portfolio

US Endoscopy has been a trusted partner in foreign body and polyp retrieval for over 20 years. Our broad portfolio of 20 different Roth Net® retrievers consists of sterile and non-sterile devices in a variety of sizes and shapes supporting polyp retrieval to the toughest foreign body and food bolus cases.

	Roth Net® retriever - standard	Roth Net® Platinum® retriever - food bolus	Roth Net® Platinum® retriever - universal	Roth Net® retriever - select	Roth Net® retriever - polyp	Roth Net® retriever - mini polyp	Roth Net® retriever - maxi	Roth Net® retriever - mini	Roth Net® retriever - enteroscope
	BX00711050 / 1052 / 1053 / 1162	BX00711155	BX00715050 / 5060	BX00711169 / 1170	BX00711150 / 1151	BX00711158 / 1159	BX00711059 / 1165	BX00711057 / 1168	BX00711152
Device									
Objects									
	FOREIGN BODY RETRIEVAL				POLYP RETRIEVAL		SPECIALTY RETRIEVAL		
	Round/blunt foreign bodies	Food bolus	Most foreign bodies	Food bolus/most foreign bodies	Polyps	Polyps	Large foreign bodies/ EMR specimens	Small foreign bodies	Foreign bodies/polyps in small bowel
Size	Net size: 3.0cm x 6.0cm Length: 230cm	Net size: 4.0cm x 5.5cm Length: 160cm	Net size: 4.0cm x 5.5cm Length: 230cm	Net size: 3.0cm x 6.0cm Length: 230cm	Net size: 3.0cm x 6.0cm Length: 230cm	Net size: 2.5cm x 4.5cm Length: 230cm	Net size: 4.0cm x 8.0cm Length: 160cm / 230cm	Net size: 2.0cm x 4.5cm Length: 160cm	Net size: 3.0cm x 6.0cm Length: 350cm
Features	<ul style="list-style-type: none"> Strong/reinforced net 	<ul style="list-style-type: none"> Large/octagonal shape Opens to maximum capacity over multiple passes 	<ul style="list-style-type: none"> Large/octagonal shape Strong/reinforced net for durability 	<ul style="list-style-type: none"> Alternative net-weave for enhanced durability Unique netting for improved visibility 	<ul style="list-style-type: none"> More transparent netting Easily encapsulates excised polyps 	<ul style="list-style-type: none"> Round shape with memory points for consistent shape after multiple passes 	<ul style="list-style-type: none"> Largest capacity retrieval net on the market Specifically designed for large foreign bodies and polyps 	<ul style="list-style-type: none"> Designed for smaller foreign bodies Compatible with smaller channel endoscopes 	<ul style="list-style-type: none"> For small bowel retrieval

Please contact us today for continuing education and in-servicing offerings. For more information on our Foreign Body Management solutions, visit www.usendoscopy.com/foreignbody.

Our products are proudly: **made in the USA**

US Endoscopy is a subsidiary of STERIS Corporation.
©2017 US Endoscopy. All Rights Reserved.
Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy

US endoscopy
listening...and delivering solutions®
www.usendoscopy.com