











cleaning brush reference chart

	Channel					Channel and Valve	Valve	Specialty		
										
	Standard Channel Cleaning Brush	Jumbo Channel Cleaning Brush	Jumbo Channel Cleaning Brush - Enteroscopy	Gemini Brush™ double-ended channel brush	Wooly Brush® channel brush	Double-Header™ combination channel and valve/control head cleaning brush	DuoSwift™ combination squeegee brush	Valve/Control Head Cleaning Brush	Small Wonder® channel cleaning brush	Small Wonder® channel cleaning brush
Product code	BX00711601	BX00711603	BX00711606	BX00711607 BX00711627	BX00711612	BX00711609 BX00711629	BX00711618 BX00711619	BX00711610	BX00711616	BX00711617
Shaft length	230cm	230cm	350cm	230cm	230cm	230cm	230cm	—	160cm	230cm
Catheter stiffness	standard	standard	standard	standard	stiff	standard	standard	—	specialty	specialty
Brush head	standard	large	large	standard	standard longer length	standard valve/control head	squeegee/standard	valve/control head	specialty	specialty
Cleaning method	push	push	push	pull	push/pull	push	push	push	push	push
Channel size	2.0-7.0mm	2.8-9.0mm	2.8-9.0mm	2.0-7.0mm	2.0-7.0mm	2.0-7.0mm	2.8-7.0mm	—	1.0-2.25mm	1.0-2.25mm

US Endoscopy is a subsidiary of STERIS Corporation.
 © US Endoscopy or its affiliates. All rights reserved.
 Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy