



Oracle[®] EUS balloons



listening...and delivering solutions[®]

Broadening our portfolio of products

US Endoscopy is proud to expand the Oracle® EUS balloon portfolio to include two new models. These balloons are designed to provide a clear view throughout endoscopic ultrasound procedures for tissue acquisition from lesions.

The Oracle® EUS non-latex balloon...

- has been manufactured not to include latex
- elastic band design allows for easy application and offers a secure grip to ensure the balloon stays on the endoscope
- is compatible with both Radial and Curvilinear echoendoscopes



Radial array



Curvilinear array

product image	product number	description	sterile	units/box
	BX00711891	Oracle® EUS non-latex balloon – Olympus curvilinear / radial array Fujifilm curvilinear*	yes	20

*This product has been cleared for use in endoscopic procedures that utilize a balloon in conjunction with an echoendoscope in the upper and lower gastrointestinal tract. The balloon is compatible with Olympus 140, 240, 160, 260 and Fujifilm curvilinear array echoendoscopes (models EG-530UT and EG-530UT2).



ect solutions for EUS procedures...

now include both latex and non-latex options. Our Oracle® EUS balloons and procedures and complement our Moray® micro forceps, which is ideal for use within and outside the GI tract.

The Oracle® EUS latex balloons...

- provide a quality image
- are cost effective options
- are available for both Radial and Curvilinear echoendoscopes
- have been manufactured to contain latex



product image	product number	package color	description	sterile	units/box
	BX00711892	navy blue	Oracle® EUS latex balloon - Olympus radial	yes	20
	BX00711893	baby blue	Oracle® EUS latex balloon - Olympus linear*	yes	20

*This product has been cleared for use in endoscopic procedures that utilize a balloon in conjunction with an echoendoscope in the upper and lower gastrointestinal tract. The balloons are compatible with Olympus 140, 240, 160, 260 and 180 curvilinear array echoendoscopes.

Please contact your US Endoscopy representative today for more information on the Oracle® EUS balloons.

Distributed by:



Olympus® is a registered trademark of Olympus Optical, Inc.
US Endoscopy is a subsidiary of STERIS Corporation.
Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks
denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.
©2017 US Endoscopy. All Rights Reserved.

5976 Heisley Road
Mentor, OH 44060
phone (1) 440 / 639-4494
fax (1) 440 / 639-4495
customer service (1) 800 / 769-8226

www.usendoscopy.com