

The Dream Team...

of foreign body management devices



The Dream Team...

For over 25 years, STERIS | US Endoscopy has led the innovation and development of the specialty tools needed for successful Foreign Body retrieval. We are proud to offer a complete and extensive portfolio of retrieval solutions to ensure you have the right product for all clinical situations.



Guardus[®] overtube

Gentle protection and enhanced safety during foreign body extraction.

BX00711146 (25cm length)
BX00711148 (50cm length)
standard scopes (8.6-10.0mm scope OD)

BX00711147 (25cm length)
BX00711149 (50cm length)
therapeutic scopes (9.8-11.7mm scope OD)

Raptor[®] grasping device

Offers the versatility needed for retrieving a variety of foreign bodies and stents.

BX00711177
sheath diameter 2.4mm
length 230cm

BX00711178
sheath diameter 1.9mm
length 200cm

Talon[®] grasping device

Grips into the food impaction for a swift retrieval.

BX00711175
sheath diameter 2.5mm
length 160cm

Roth Net® retrievers

Our broad portfolio of Roth Net® retrievers supports a variety of clinical applications and retrieval needs. From foreign bodies to challenging food bolus cases, the following will cover you from top to bottom.

NEW Roth Net® retriever - 360

BX00711197

net size 3.0cm x 5.0cm
sheath diameter 2.5mm
length 230cm

Opens to full capacity after multiple passes




Roth Net® Platinum® retriever - universal

BX00715050

net size 4.0cm x 5.5cm
sheath diameter 2.5mm
length 230cm

Opens to full capacity after multiple passes



Roth Net® retriever - select

BX00711170

net size 3.0cm x 6.0cm
sheath diameter 2.5mm
length 230cm



Roth Net® retriever - maxi

BX00711059

net size 4.0cm x 8.0cm
sheath diameter 3.0mm
length 160cm



Roth Net® retriever - mini

BX00711057 / BX00711168*

net size 2.0cm x 4.5cm
sheath diameter 1.8mm
length 160cm



Sterile options available.



Falcon® rotatable retrieval basket

A great utility tool to have on hand for removal of long objects and other foreign bodies not suitable for the Roth Net® retriever.

BX00711141

basket size 3.5cm x 6.0cm
sheath diameter 2.5mm
length 230cm

STERIS | US Endoscopy offers Continuing Education and in-servicing to ensure you and your staff are prepared for any clinical challenge – contact us today to schedule.

www.usendoscopy.com/education



5976 Heisley Road
Mentor, OH 44060
phone (1) 440 / 639-4494
fax (1) 440 / 639-4495
customer service (1) 800 / 769-8226
email info@usendoscopy.com

www.usendoscopy.com

US Endoscopy is a subsidiary of STERIS Corporation.

©2019 US Endoscopy. All Rights Reserved.

All marks denoted with © or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.