



listening...and delivering solutions®

5976 Heisley Road
Mentor OH 44060

phone 440 / 639.4494

fax 440 / 639.4495

customer service 800 / 769.8226

www.usendoscopy.com

AdvanCE® capsule endoscope delivery device Performance and Compatibility Guarantee

US Endoscopy stands behind the performance and compatibility of its products.

All US Endoscopy products, including the AdvanCE® device, are designed and manufactured in accordance with the FDA's Quality Systems Regulations and International Standards (ISO 13485). Extensive testing performed by US Endoscopy has shown the AdvanCE® device to be compatible with the vast majority of video capsules manufactured by Given Imaging, Olympus, and IntroMedic.

Please refer to the following table for video capsule compatibility:

Product Description	Product Number	Sheath Diameter	Sheath Length	Compatible Capsule Dimensions
AdvanCE® capsule endoscope delivery device	BX00711144	2.5mm	180cm	Diameter: 10.5mm-11.5mm Length: 23.5mm-26.5mm

If a customer is not satisfied with the AdvanCE® device for any reason, US Endoscopy offers a 100% refund on all unused and unopened product for one year from the date of shipment. Please refer to US Endoscopy's terms and conditions for additional details.

The AdvanCE® capsule endoscope delivery device allows for direct placement of video capsules into the stomach or duodenum in patients who are unable to swallow or pass the capsule through the pylorus. Controlled endoscopic placement using the AdvanCE® capsule endoscope delivery device helps minimize the risk of stalled capsules in the esophagus or stomach related to anatomical or motility disorders.

The AdvanCE® device is to be used by trained personnel who are familiar with the instructions for use of the product. Please contact your US Endoscopy Product Specialist if you have any questions on the use of this product.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or TM are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.