

AdvanCE[®] capsule endoscope delivery device



Direct placement of video capsules...

The *AdvanCE® capsule endoscope delivery device* allows for direct endoscopic placement of video capsules in patients who are unable to swallow or pass the capsule through the pylorus.



The AdvanCE® device...

- Facilitates endoscopic placement of video capsules into the stomach or duodenum
- Minimizes risk of stalled capsule in esophagus or stomach, related to anatomical or motility disorders.



For more information on the AdvanCE® delivery device, contact US Endoscopy today.

product number	sheath diameter (mm)	length (cm)	compatible video capsule dimensions (mm)	units/box
BX00711144	2.5	180	diameter: 10.5 - 11.5 length: 23.5 - 26.5	3

Distributed by:



listening...and delivering solutions®

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639.4494

fax +1 440 / 639.4495

email global@usendoscopy.com

www.usendoscopy.com

760398 Rev. C

US Endoscopy is a subsidiary of STERIS Corporation.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

©2017 US Endoscopy or its affiliates. All rights reserved.

Follow us:    