

# Take-Out® transport pad



# A safe and easy way to transport your endoscope...

The Take-Out® transport pad is designed for easy setup and to completely enclose and transport endoscopes, instruments, and accessories to the reprocessing room after use.

## The Take-Out® transport pad offers...



- Simple setup of equipment and procedural supplies
- Absorbent pad helps prevent leaking fluids and protect work area
- Transparency allows for visualization of endoscope
- Drawstrings completely **enclose** soiled endoscope within pad for safe transport
- Enclosed pad provides protection for endoscope and staff
- Included biohazard sticker allows for proper labeling

For more information on the Take-Out® transport pad, contact US Endoscopy today.

| product number | description   | dimensions  | units/ box |
|----------------|---------------|-------------|------------|
| BX00711760     | transport pad | 98cm x 75cm | 25         |

Distributed by:



5976 Heisley Road  
Mentor, OH 44060 USA

phone +1 440 / 639,4494

fax +1 440 / 639,4495

email [global@usendoscopy.com](mailto:global@usendoscopy.com)

[www.usendoscopy.com](http://www.usendoscopy.com)

760842 Rev. C

US Endoscopy is a subsidiary of STERIS Corporation.  
All marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.  
©2018 US Endoscopy. All Rights Reserved.

Follow us:    