

AquaShield® water bottle system




listening...and delivering solutions®

Proactively support patient safety...

The disposable AquaShield® system offers an alternative option to reprocessing reusable water bottle systems.



AquaShield® system consists of cap, tubing, connector, & 24-hour use sticker



The disposable AquaShield® system...

- Offers a convenient solution, saving time and money associated with reusable systems
- Reduces the risk of cross-contamination between patients
- A 24-hour use system - use one for the day and throw it away
- Disposable solution for Olympus, Pentax and Fujifilm users

For more information on the AquaShield® water bottle system, contact US Endoscopy today.

| product number | description | CO ₂ compatibility | endoscope compatibility | units/box |
|----------------|--|-------------------------------|-------------------------|-----------|
| BX00711540 | AquaShield® system | no | Olympus 140/160/180/190 | 25 |
| BX00711541 | AquaShield® system - CO ₂ (4" line) | yes | Olympus 140/160/180/190 | 25 |
| BX00711549 | AquaShield® system - CO ₂ (60" line) | yes | Olympus 140/160/180/190 | 25 |
| BX00711550 | AquaShield® system - Pentax | yes | Pentax | 25 |
| BX00711551 | AquaShield® system CO ₂ - Pentax (60" line) | yes | Pentax | 25 |
| BX00711553 | AquaShield® system CO ₂ - Fujifilm (60" line) | yes | Fujifilm | 25 |

NEW
NEW
NEW

All AquaShield® systems are compatible with commercially available bottles such as Baxter/Hospira/B.Braun water bottles. AquaShield® system CO₂ caps are capable of being used with air or CO₂ insufflation.

B. Braun refers to B. Braun Medical Inc.
Baxter refers to Baxter International, Inc.
Hospira refers to Hospira, Inc.
Olympus® is a registered trademark of Olympus Corporation.
Fujifilm® is a registered trademark of Fujifilm Corporation.
PENTAX® is a registered trademark of Hoya Corporation.
US Endoscopy is a subsidiary of STERIS Corporation.
©2017 US Endoscopy. All Rights Reserved.

Listening...and delivering solutions® the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.



listening...and delivering solutions®

5976 Heisley Road
Mentor, OH 44060

phone 440 / 639.4494

fax 440 / 639.4495

customer service 800 / 769.8226

www.usendoscopy.com

Follow us:

