

AdvanCE[®] capsule endoscope delivery device



listening...and delivering solutions[®]

Direct placement of video capsules...

The *AdvanCE*[®] capsule endoscope delivery device

allows for direct endoscopic placement of video capsules in patients who are unable to swallow or pass the capsule through the pylorus.



The *AdvanCE*[®] system...

- Facilitates endoscopic placement of video capsules into the stomach or duodenum
- Minimizes risk of stalled capsule in esophagus or stomach, related to anatomical or motility disorders.

For more information on the AdvanCE[®] delivery system, contact US Endoscopy today.

product number	sheath diameter (mm)	length (cm)	compatible video capsule dimensions (mm)	units/box
BX00711144	2.5	180	diameter: 10.5 - 11.5 length: 23.5 - 26.5	3



listening...and delivering solutions[®]

5976 Heisley Road
Mentor OH 44060

phone 440 / 639.4494

fax 440 / 639.4495

customer service 800 / 769.8226

www.usendoscopy.com

US Endoscopy is a subsidiary of STERIS Corporation.

Listening...and delivering solutions[®], the US Endoscopy listening...and delivering solutions design, and all marks denoted with © or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

©2017 US Endoscopy or its affiliates. All rights reserved.

Follow us:    