




**ArC Smart™** argon coagulation probe  
& ArConnect® argon probe connector



*listening...and delivering solutions®*

## The perfect companion to the superior performance of the gi4000 ArC Smart™ argon linear beam...

The ArC Smart™ argon coagulation probe ensures optimum beam and energy delivery when needed.

### The ArC Smart™ argon coagulation probe offers...

- Quick ignition and a stable beam
- Ideal stiffness for flexible endoscopic maneuvers, including side viewing scopes
- Availability in multiple sizes and lengths



*One single-use, ArConnect® probe connector is required for each argon coagulation procedure in order to attach standard connect compatible argon probes.*

For more information on the ArC Smart™ argon coagulation probe and specialty consumables, contact US Endoscopy today.

product number	description	sheath diameter	length	units/box
G1140001	ArC Smart™ probe - straight fire	2.3mm	220cm	10
G1140002	ArC Smart™ probe - straight fire	2.3mm	330cm	10
G1140003	ArC Smart™ probe - straight fire	3.2mm	220cm	10
G1140004	ArC Smart™ probe - straight fire AND ArConnect® probe connector	2.3mm	220cm	10 each / 20 pieces
G1140005	ArC Smart™ probe - straight fire AND ArConnect® probe connector	2.3mm	330cm	10 each / 20 pieces
G1140006	ArC Smart™ probe - straight fire AND ArConnect® probe connector	3.2mm	220cm	10 each / 20 pieces
G1120020	ArConnect® probe connector	/	/	10



*listening...and delivering solutions®*

5976 Heisley Road  
Mentor, OH 44060

phone 440 / 639.4494

fax 440 / 639.4495

customer service 800 / 769.8226

[www.usendoscopy.com](http://www.usendoscopy.com)

US Endoscopy is a subsidiary of STERIS Corporation.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

©2017 US Endoscopy. All Rights Reserved.

Follow us:    