



**Gonogo™**  
monopolar active cord tester



*listening...and delivering solutions®*

## Is your procedure a go?

The Gonogo™ monopolar active cord tester allows you to trouble-shoot your active cords to check for functionality.

### The Gonogo™ monopolar active cord tester...

- features a simple design – insert both ends of the active cord into the test ports and push the black button. Green light for Go and red light for No Go.
- offers compatibility with most leading endoscopy active cords.
- provides simple, reliable testing - the active cord does not need to be connected to the electro-surgical unit to test.



For more information on the Gonogo™ monopolar active cord tester, contact US Endoscopy today.

product number	description	units/box
G1609	Gonogo™ monopolar active cord tester	1
G11-120-06	Active cord (monopolar, 3mm) 10ft, narrow handle, blue	1

US Endoscopy is a subsidiary of STERIS Corporation.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

©2017 US Endoscopy. All Rights Reserved.

Follow us:



**us endoscopy**

*listening...and delivering solutions®*

5976 Heisley Road  
Mentor, OH 44060

phone 440 / 639.4494

fax 440 / 639.4495

customer service 800 / 769.8226

[www.usendoscopy.com](http://www.usendoscopy.com)

761427 Rev. A