

# Disposable BioShield™ Biopsy Valves & Irrigation Accessories

Patent Pending

US Endoscopy, a world leader in endoscopy device design and manufacturing, introduces a new line of Disposable BioShield™ Biopsy Valves and Irrigation Accessories to minimize leakage and back-flow of fluid and other biomaterial from the endoscope's accessory channel. The design also maximizes intraprocedural irrigation capability. All accessories are disposable and offer unsurpassed performance, value and protection to caregivers from unnecessary exposure to biomaterial.

## Universal Irrigating Adaptor

Atraumatic, blunt end adaptor attaches to luer-lock and slip-tip syringes. Maximizes direct, intraprocedural irrigation. Compatible with a wide variety of endoscope biopsy valves.

## BioShield Biopsy Valve with Universal Irrigating Adaptor

Engineered to minimize biomaterial leakage from the biopsy port throughout the procedure. Accommodates a variety of different sized instruments. Press fit design incorporates a reinforced, thicker wall allowing effortless instrument passage through the biopsy port while maintaining structural integrity. Includes Universal Irrigating Adaptor.

*Recommended for upper GI procedures*

## BioShield Biopsy Valve with Irrigation Line

Same features as the Biopsy Valve. Irrigation line allows remote access for irrigation. Proximal one-way valve prevents back-flow of irrigation fluids.

*Recommended for lower GI procedures*

**711131**  
Universal  
Irrigating  
Adaptor  
*(syringe not  
included)*

**711132**  
Biopsy Valve  
w/Universal  
Irrigating  
Adaptor

**711133**  
Biopsy Valve  
w/Irrigation  
Line  
*(syringe not  
included)*

### US Endoscopy Disposable BioShield Biopsy Valves & Irrigation Accessory Specifications

Product Number	Description	Adaptor Gauge	Adaptor Length	Irrigation Line Length	Units/Box
711131	Universal Irrigating Adaptor	15	1"	-	50
711132	Biopsy Valve w/Universal Irrigating Adaptor – Olympus® scopes <i>Recommended for upper GI procedures</i>	15	1"	-	50
711133	Biopsy Valve w/Irrigation Line – Olympus® scopes <i>Recommended for lower GI procedures</i>	-	-	12"	50

Olympus is a registered trademark of Olympus Optical Co., Ltd.

US Endoscopy biopsy valve & irrigation accessories are intended to be used with gastrointestinal endoscopes.



*listening...and delivering solutions*

5976 Heisley Road  
Mentor OH 44060

phone 440 / 639.4494

fax 440 / 639.4495

customer service 800 / 769.8226

[www.usendoscopy.com](http://www.usendoscopy.com)

**Ordering Information 800/769-8226**