

“Pull” PEG Systems



US Endoscopy’s initial placement 20FR and 24FR “Pull” silicone feeding catheter systems are the most up-to-date devices on the market today. These devices are unparalleled in quality, and were specifically designed with the help of doctors and nurses to promote efficient, hassle-free PEG and PEJ placements. The US Endoscopy Pull PEG System combines procedure tray and support tray items, as well as a deluxe grasping snare.



- **Easy traction removal**
Can be removed without endoscopic assistance
- **Deluxe grasping snare included**
Convenient; addresses cost containment
- **Centimeter marked catheter**
Helps monitor insertion depth and migration
- **Soft, flexible silicone dome**
Protects against leakage and stoma irritation
- **Nylon coated “pull” wire**
Minimizes friction
- **One piece silicone external bolster**
Speeds bolster positioning
- **20FR and 24FR sizes available**
Optimal for gastric feeding and compatible with 9FR and 12FR J-Tubes
- **Instructions for Use, Patient Labels and Patient Care Instruction Manual attached to outside of each tray**
Easy accessibility to the items you need first



5976 Heisley Road, Mentor, Ohio 44060
Phone (440) 639-4494 • Fax (440) 639-4495
Customer Service (800) 769-8226
Email: info@usendoscopy.com www.usendoscopy.com



ORDERING INFORMATION 800-769-8226

Product Number	Description	Feeding Tube Length	Feeding Tube Size
00710801	Pull PEG System	60cm	20 FR
00710805	Pull PEG System	76cm	24 FR

Includes:

- 1 20FR 60 cm or 24FR 76 cm Silicone Gastrostomy Feeding Tube, with insertion depth markings and pull wire
- 1 Silicone External Bolster, One Piece
- 1 Silicone Dual Port Feeding Adaptor
- 1 Nylon Coated Pull Wire
- 1 Grasping Snare
- 1 24" x 36" Fenestrated Drape
- 1 No. 11 Scalpel
- 1 5cc Syringe
- 1 Filter Cannula
- 1 25 Gauge x 5/8" Needle
- 1 19 Gauge x 1 1/2" Needle
- 1 14 Gauge Introducer Catheter
- 1 5cc Ampule Lidocaine
- 1 Triple Pack Povidone-Iodine Swabs
- 1 Prep Razor
- 4 4" x 4" Gauze Pads
- 4 2" x 2" Cross Cut Drainage Gauze
- 2 Ointment (10% Povidone-Iodine)
- 2 Lubricating Jelly
- 1 Curved Hemostat
- 1 Stainless Steel Scissors
- 2 Cotton Tip Applicators
- 1 48FR Bite Block
- 1 Easy Access Patient Care Manual
- 1 Easy Access Triple Patient Label
- 1 Easy Access Instructions For Use

