


# Disposable Endoscope Bite Blocks



US Endoscopy offers an expansive line of bite blocks to meet the needs of a variety of patients. **All materials used to produce our bite blocks are LATEX FREE** in order to help you safely serve the growing population of latex sensitive patients and healthcare staff. US Endoscopy bite blocks can be purchased with or without a dental retention rim\*, and are available in 38FR (infants and small children), 48FR (school age children through geriatrics), and 60FR (adult) sizes. All sizes offer lateral side ports for suctioning and finger-assisted scope insertion. US Endoscopy bite blocks are made of a sturdy yet forgiving material that will protect your delicate instruments, but offer comfort for your patients. Their flexible frame molds gently to your patients' faces - especially nice for babies and toddlers with chubby cheeks!

\* The dental retention rim helps keep bite blocks in place in healthy adolescents and adults, but can be cumbersome for infants, toddlers, geriatrics and patients with limitations of the mouth or jaw to accommodate. US Endoscopy offers bite blocks both with and without a dental retention rim to serve all of your patients' needs.



- **All materials and retaining straps are LATEX FREE**  
*Safely serve the growing population of latex sensitive patients and staff*
- **Offered with a Latex free rubber strap, Velcro Strap or without a retaining strap**
- **Available with or without a dental retention rim**  
*Addresses a variety of patients, including those with limitations of the mouth or jaw*
- **Multiple sizes and styles**  
*Addresses patients of all ages and clinician preference*
- **Realistically priced**  
*No need to worry about extra cost associated with Latex free upgrade*



5976 Heisley Road, Mentor, Ohio 44060  
Phone (440) 639-4494 • Fax (440) 639-4495  
Customer Service (800) 769-8226  
Email: [info@usendoscopy.com](mailto:info@usendoscopy.com) [www.usendoscopy.com](http://www.usendoscopy.com)



## ORDERING INFORMATION 800-769-8226

Product Number	Description	Opening Size	Side Port	Dental Retention
00712801	Standard Bite Block - Latex Free Strap	16mm (48FR)	No	Yes
00712802	Standard Bite Block - No Strap	16mm (48FR)	No	Yes
00712803	Standard Bite Block - Velcro Strap	16mm (48FR)	No	Yes
00712804	Maxi Bite Block - Latex Free Strap	20mm (60FR)	Yes	Yes
00712805	Maxi Bite Block - No Strap	20mm (60FR)	Yes	Yes
00712806	Maxi Bite Block - Velcro Strap	20mm (60FR)	Yes	Yes
New! 00712808	Pediatric Bite Block - Latex Free Strap	12.6mm (38FR)	Yes	No
New! 00112809	Intermediate Bite Block - Latex Free Strap	16mm (48FR)	Yes	No

