


US Endoscopy flexible Ponsky Cholangiocatheters are 5FR catheters that offer the option of a .025 inch coated guidewire for **easier cystic duct cannulation**. These catheters have bold markings on their distal tips to identify insertion depth, and are compatible with cholangiocatheter clamps. For those who prefer not to use a guidewire, the catheters are also available with a **round or taper tip**.



- **Two basic catheter tip configurations; Taper, Round tip**
Surgeon preference; select the tip that best fits the anatomical findings
- **Guidewire**
Allows for easier cannulation of the cystic duct
- **Clear catheter**
Allows surgeon to identify the presence of air bubbles when contrast media is injected
- **Black insertion markings**
Aids in determining depth of catheter insertion
- **Realistically priced**
Addresses cost containment issues

ORDERING INFORMATION 800-769-8226

Product Number	Catheter Diameter	Catheter Length (cm)	Distal Tip Configuration	Distal Tip Diameter
00712031	5 FR	68	Round w/ Guidewire	5 FR
00712032	5 FR	68	Taper	4 FR

