

# Articulator™ injection needles



**US** endoscopy

*listening... and delivering solutions®*

# Injection therapy as it should be...

The Articulator™ injection needle from US Endoscopy™ is a reliable, yet affordable needle, ideal for everyday use and dependable enough for your most immediate needs. From saline assisted polypectomy, to esophageal varices, ulcers, and tattooing, the Articulator™ needle is the steadfast choice for injection therapy.

## *The Articulator™ needle offers a...*

- *Luer-lock, spring-loaded handle to determine if the needle is either engaged or retracted, no in-between – reducing the risk of inadvertent scope damage and incidental needle sticks, while providing maximum needle projection*
- *Kink resistant, stainless steel sheath to pass smoothly and perform in fully articulated endoscopes*
- *Stainless steel hub to aid in tamponade and help ensure only the needle penetrates the mucosa*

For more information on Articulator™ needles, contact US Endoscopy today.

product number	sheath diameter (mm)	length (cm)	needle projection (mm)	needle gauge	units/box
00711807	2.5	230	5	25	5
00711808*	2.5	350	5	25	5

\* 00711808 does not have the spring-loaded handle.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

**us endoscopy**

*listening...and delivering solutions\**

5976 Heisley Road  
Mentor OH 44060

phone 440 / 639.4494

fax 440 / 639.4495

customer service 800 / 769.8226

[www.usendoscopy.com](http://www.usendoscopy.com)

760422 Rev. A