

bite blocks

disposable and latex-free



us endoscopy

listening... and delivering solutions®

An impressive selection – latex-free and worry free

The US Endoscopy family of bite blocks includes a variety of sizes, designs, and models to ensure the best fit for the patient and clinician. Sturdy, yet forgiving material protects delicate instruments and offers comfort for your patient.

Our bite blocks and coordinating straps are also disposable and latex-free, addressing latex sensitivity and risk of severe allergic reaction. 



selections with a dental retention rim

- Dental retention rim helps to keep bite block in place

selections without a dental retention rim

- Designed with a smooth insertion rim and flexible frame to mold gently to your patient's face, these latex free bite blocks are especially suited for those with limitations of the mouth or jaw



Call US Endoscopy today to learn more about our family of bite blocks

product number	description	opening size	units/box
with dental retention rim			
00712801	standard – latex-free strap	16mm / 48 FR	100
00712802	standard – no strap	16mm / 48 FR	100
00712803	standard – velcro strap	16mm / 48 FR	100
00712804	maxi bite – latex-free strap	20mm / 60 FR	100
00712805	maxi bite – no strap	20mm / 60 FR	100
00712806	maxi bite – velcro strap	20mm / 60 FR	100
without dental retention rim			
00712808	pediatric – latex-free strap	12,6mm / 38 FR	25
00712809	intermediate – latex-free strap	16mm / 48 FR	100

us endoscopy

*listening...and delivering solutions**

5976 Heisley Road
Mentor OH 44060

phone 440 / 639.4494

fax 440 / 639.4495

customer service 800 / 769.8226

www.usendoscopy.com