


eTrap[®] polyp trap



listening...and delivering solutions[®]

Inspired by nurses...

When US Endoscopy set out to design the ideal polyp trap, we interviewed numerous clinicians, formulating an impressive wish list. The result: a polyp trap that allows for quick and accurate retrieval, while safeguarding clinicians from unnecessary exposure to biomaterial. The eTrap[®] polyp trap offers the simplicity desired by clinicians and effectiveness both you and your patients deserve.

The eTrap[®] - a polyp trap designed with you in mind...

Removable specimen strainers provide:

- Little or no biomaterial exposure - no more sifting through fluid to locate a specimen
- Simple transport and preservation of intact polyp specimens - dip the tray in the specimen jar or brush with your fingertip, and the specimen will slide right out
- Ability to collect multiple polyps with uninterrupted suction - just alternate the trays



Collection chamber features:

- Magnifying window to easily visualize the collected specimen
- Clean design that cycles fluid into the main suction receptacle



For more information on the eTrap[®] polyp trap, contact US Endoscopy today.

product number	description	units/box
00711099	polyp trap with 2 specimen trays/unit	25 units/box

us endoscopy

listening...and delivering solutions[®]

5976 Heisley Road
Mentor OH 44060

phone 440 / 639.4494

fax 440 / 639.4495

customer service 800 / 769.8226

www.usendoscopy.com

Listening...and delivering solutions[®], the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy

760264 Rev B