

Velocity[®]
irrigation system




Experience the power of irrigation...

Improve your irrigation capabilities through the *Velocity® irrigation system*.

The system provides critical elements required for optimal irrigation:

- Increases irrigation volume
- Provides variable pressure
- Allows for immediate physician control



Intraprocedural irrigation using Velocity® system

The system is an ideal solution for cases that require ample irrigation capabilities:

- Poor preparations in the colon
- Gastrointestinal bleeds
- Water colon insufflation
- Blood clot removal
- Undigested food in the stomach

The *Velocity® irrigation system* includes an easy-to-use manual foot pump and irrigation tubing that connects directly to your endoscope and sterile water supply.



Easy-to-use manual foot pump



Velocity® system connects directly to your endoscope



Irrigation tubing connects to sterile water supply (water bottle not included)

For more information on the Velocity® irrigation system, contact US Endoscopy today.

product number	description	scope compatibility	sterile	units/box
00711520**	Velocity® irrigation pump	n/a	yes	5
00711521	Velocity® irrigation pump	n/a	yes	5
00711139	Velocity® biopsy valve	Olympus/Fujinon*	no	25
00711140	Velocity® biopsy valve	Pentax	no	25

Distributed by:



*listening...and delivering solutions**

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639,4494

fax +1 440 / 639,4495

email global@usendoscopy.com

www.usendoscopy.com

760346 Rev. E

*G5 series or newer

**Confirm product number with the distributor of US Endoscopy's products

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy

Fujinon® is a registered trademark of Fuji Film Co., Ltd.

Olympus® is a registered trademark of Olympus Optical, Inc.

PENTAX® is a registered trademark of Asahi Kogaku Kogyo Kabushiki Kaisha