

Carr-Locke injection needle



US endoscopy

*listening... and delivering solutions**

Trusted performance...

World renowned for its reliability and consistency, the *Carr-Locke injection needle* is ideal for your everyday procedures, as well as your most tortuous, difficult injection therapy maneuvers, including use in duodenoscopes.

now available in 23 gauge

- *Stainless steel, Teflon coated, spring sheath* to virtually eliminate sheath kinking, while offering exceptional functionality and needle performance during challenging procedures
- *Ultra smooth distal metal hub* aids in tamponade and is gentle on friable tissue, helping to prevent sheath perforations
- *Luer-lock, spring-loaded handle* automatically retracts the needle, reducing the risk of inadvertent needle sticks or scope damage

Consistently deploys in retroflexed position

For more information on the Carr-Locke injection needle, contact US Endoscopy today.

product number	sheath diameter (mm)	length (cm)	needle projection (mm)	needle gauge	units/box
00711811	2.5	230	5	25	5
00711812	2.5	230	5	23	5
00711813 (Pentax)	2.5	230	5	25	5

New

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy
Teflon® is a registered trademark of Dupont

US endoscopy

listening...and delivering solutions®

5976 Heisley Road
Mentor, OH 44060

phone 440 / 639.4494

fax 440 / 639.4495

customer service 800 / 769.8226

www.usendoscopy.com

760360 Rev. C