









# essential foreign body equipment

	Roth Net® Platinum® food bolus retrieval net	Roth Net® maxi retrieval net	Roth Net® pediatric retrieval net	Talon™ grasping device	Raptor™ grasping device	Guardus® esophageal/gastric overtube	rotatable retrieval basket	Rotator™ rotatable snare
								
	0071 <b>1155</b> Net size - 4.0cm x 5.5cm	0071 <b>1059</b> Net size - 4.0cm x 8.0cm	0071 <b>1057</b> Net size - 2.0cm x 4.5cm	0071 <b>1175</b> Length - 160cm	0071 <b>1177</b> Length - 230cm	0071 <b>1146/0071 1147</b> 0071 <b>1148/0071 1149</b> Length - 25cm/50cm	0071 <b>1141</b> Basket size 3.5cm x 6.0cm	0071 <b>1120</b> Snare size - 2.5cm x 54cm
Larger blunt object removal (AA battery..)		✓				✓		
Smaller blunt object removal (buttons, coins..)	✓	✓	✓		✓	✓		
Long object removal (spoon, toothbrush..)						✓	✓	✓
Sharp-pointed object removal (toothpicks, bones, paperclips..)				✓	✓	✓	✓	
Disk battery removal	✓	✓	✓			✓		
Food bolus removal	✓			✓		✓		

## are you prepared?

**The ASGE Guideline** for the management of ingested foreign bodies states, "Equipment that should be readily available includes **rat tooth and alligator forceps, polypectomy snare, polyp grasper, Dormier basket, retrieval net, overtubes of esophageal and gastric lengths, and a foreign body protector hood.**"

*Gastrointestinal Endoscopy 11/2002*

*The above chart suggests applications for US Endoscopy foreign body management devices. Only persons having adequate training and proficiency in endoscopic foreign body removal technique are qualified to determine appropriate device usage in any given procedure.*

product	qty.	complete foreign body inventory checklist														
Roth Net® retrieval nets																
Talon™ grasping device																
Guardus® esophageal overtubes																
Guardus® gastric overtubes																
Rotator™ rotatable snare																
rotatable retrieval basket																
Raptor™ grasping device																
polyp grasper																
foreign body protector hood																
<b>initials &amp; date</b>																

product	qty.	complete foreign body inventory checklist														
Roth Net® retrieval nets																
Talon™ grasping device																
Guardus® esophageal overtubes																
Guardus® gastric overtubes																
Rotator™ rotatable snare																
rotatable retrieval basket																
Raptor™ grasping device																
polyp grasper																
foreign body protector hood																
<b>initials &amp; date</b>																

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.



phone (1) 440 / 639-4494  
 fax (1) 440 / 639-4495  
 customer service (1) 800 / 769-8226  
[www.usendoscopy.com](http://www.usendoscopy.com)